

Why should you use the Activating Cleansers with your LumiSpa device?

There are five ageLOC LumiSpa® Activating Cleanser options for ageLOC LumiSpa® device, each designed for different skin types: dry, normal/combo, oily, blemish prone and sensitive. This means you can tailor your LumiSpa® experience to fit your skincare regimen and have it work optimally for you.

The Activating Cleansers are formulated to achieve a certain rheology and viscosity to enable precise cushioning, comfortable grip, gliding and interaction with the skin, thus ensuring the optimal effectiveness of LumiSpa's proprietary Silicone Head. Using a different cleanser that doesn't have the same viscosity or texture as the ageLOC LumiSpa® Activating Cleansers could let the device latch onto your skin too tightly or conversely, let it slide on the face, without any additional benefits.



The Activating Cleansers are powerful products formulated using Nu Skin's proprietary ageLOC® ingredients blend and designed to work together with the patented Two-Sense Motion Technology of LumiSpa's Silicone Heads. They ensure that you receive visible improvements and cleansing benefits of the experience revealing healthy, bright and youthful-looking skin.

LumiSpa System delivers seven visible skin benefits, including softness, smoothness, radiance, clarity, freshness, reduced appearance of pores, and thoroughly clean skin. All were proven in a 12-week in vivo study, in which 150 participants were involved and noticed improvements over the course of 12 weeks using the full system (LumiSpa device + ActivatingCleanser) twice a day.

Additionally, our latest in-vivo study with 25 participants has shown that this incredible skin cleansing system also does not disturb the skin's natural microbiome. A third-party clinical laboratory was commissioned to evaluate the impact of one cleansing with five different LumiSpa combinations on the skin's natural microbiome.. Four hours post-cleansing, there was no significant change in the dominant microorganism population compared to baseline, confirming that the ageLOC LumiSpa system doesn't disturb the skin microbiome.*

How to use the **Activating** **Cleansers** with your LumiSpa device

For optimal results, use the LumiSpa Cleansing System for two minutes, both morning and evening.

In the evening, if you are wearing heavy make-up, first remove it using the Nutricentials Micellar Beauty Water or the Nu Colour Make-up Remover. If you are only wearing medium coverage to light make-up, or a sun protection, you can skip this first step and directly start using your LumiSpa. The System will efficiently remove make-up and sunscreen.



First moisten your face with lukewarm water, then place an almond-size amount of Activating Cleanser that best suits your skin type onto fingertips and apply it all over your face. Wet the silicone head and turn on your device. Glide it slowly across your face in broad strokes. Don't scrub or press too hard. The device will pause if you are misusing it.

It will pause every 30 seconds. This is your cue to move to another area of your face. After 2 minutes, your device will switch off.

Rinse your face with water and gently pat it dry, then apply any further skincare.

Applying the Activating Cleanser directly onto the LumiSpa Head can hinder its proper rotation during cleansing and cause head blockage while its drying out. Although the LumiSpa Head is made of silicon, which is easy to clean, blocked traces of Activating Cleanser in its structure can cause quicker degradation of your LumiSpa Silicone Head.

*Wondering which Activating Cleanser you should use?
Visit https://quiz.lumispacecare.com/en_NL*