



LumiSpa®
 ACCENT

Stress, age, and a busy lifestyle can take their toll—but don't let your skin give you away! ageLOC LumiSpa pairs a one-of-a-kind device that uses patented technology with ageLOC LumiSpa products. The result? Powerful treatments that help diminish the signs of aging to keep you looking radiant.



MEET LUMISPA ACCENT

The perfect accent to brighten every day. Late nights and early mornings. Stress, age, a busy lifestyle, and sun exposure. All these things can wreak havoc on the skin around your eyes. But now you can save the bags for shopping and the puffiness for warm winter jackets. Get bright, ready-for-anything eyes with ageLOC LumiSpa Accent.

The soft treatment tip gently exfoliates the skin while deeply stimulating the delicate area around the eyes. Paired with ageLOC LumiSpa IdealEyes®, Accent delivers an invigorating treatment and an instant burst of hydration, visible firming, and the wide-awake radiance you want every day. Target undereye bags, puffiness, and dark circles twice daily for a refreshed, well-rested look. All eyes on you, baby!

WHY YOU'LL LOVE IT

- Gently stimulates skin to help maximize ageLOC LumiSpa IdealEyes benefits.
- Uses proprietary LumiSpa technology to deliver ingredients to the delicate eye area while treating and stimulating skin.
- Promotes the appearance of stronger, healthier skin.
- Refines fine lines and wrinkles for a youthful look.
- The spot treatment tip is made of the softest silicone available for gentle use on delicate skin around the eyes.
- Embedded with antimicrobial silver in the silicone tip to make cleaning as easy as possible.

MEET AGELOC LUMISPA IDEALEYES

Goodbye bags, hello bright eyes. ageLOC LumiSpa IdealEyes Activating Eye Cream is specially designed to work with ageLOC LumiSpa Accent to deliver powerful benefits to the delicate skin around your eyes. Just 30 seconds on each eye

area provides immediate hydration, visible firming, brightening, and smoothing to stop the look of tired, stressed eyes in its tracks. It's the perfect routine to help with those less-than-ideal signs of a busy life around your eye area.

WHAT POWERS IT

- **Peptide Complex**—helps visibly firm and tone delicate skin.
- **Tetrahexyldecyl Ascorbate**—a stable, powerful form of vitamin C that has well-known anti-aging benefits for lines, wrinkles, and discoloration. Protects against free radical damage, brightening the eye area.
- **Plankton Extract**—helps decrease the appearance of pigmentation in the undereye area and visibly increase skin tonicity. Supports the skin's natural recovery process.
- **Squalane**—prevents moisture loss while restoring skin's suppleness and flexibility.
- **Antioxidant Complex**—protects skin against daily environmental stress.
- **ageLOC Proprietary Blend**—an anti-aging complex that targets the sources of aging to preserve the look of youth.

CLINICAL RESULTS THAT PROVE IT

Nu Skin commissioned a third-party contract research organization to evaluate the visible short- and long-term benefits of the ageLOC LumiSpa Accent system.

Clinical Evaluation

Week 2—significant improvements in the appearance of bags and puffiness.

Week 4—significant visible lifting improvements and reduction of crepey eyelids.

Week 8—significant improvements in the appearance undereye circles.

ageLOC® LumiSpa® ACCENT

Week 12—a clinical grader recorded a 98% increase in hydration in subject's skin.

HOW TO USE IT

1. Apply ageLOC LumiSpa IdealEyes® to clean skin around the eyes (not on the eyelids or eyes themselves). Turn on the device. Following the bone structure around the eye, glide LumiSpa Accent around the eye area using slow, gentle motions. Avoid placing the Accent tip directly on eyelids. Do not scrub or move across the skin quickly. Use morning and night.
2. The device will pause briefly after 30 seconds, indicating that it's time to move to the next eye area. Treat both eyes, then press the On/Off button to turn off the device.

Note: Do not use Accent directly on bare skin. Use only with ageLOC LumiSpa IdealEyes.

TRY THESE WITH IT

- ageLOC LumiSpa + Treatment Cleansers
- Tru Face® Line Corrector
- ageLOC Tru Face Essence Ultra
- ageLOC Eye Mask

LEARN MORE ABOUT IT

How does ageLOC LumiSpa Accent work?

With its soft silicone treatment tip, ageLOC LumiSpa Accent extends LumiSpa technology to move at the precise frequency to promote healthy, brighter skin around the eyes—in just one minute. The tip gently stimulates the skin around the eye, helping ageLOC LumiSpa IdealEyes nourish and treat this delicate area.

What makes ageLOC LumiSpa Accent unique?

ageLOC LumiSpa Accent is different in three key ways:

1. ageLOC LumiSpa Accent targets the delicate eye area (a top skin concern area for many consumers).
2. ageLOC LumiSpa Accent uses proprietary LumiSpa technology to gently stimulate the delicate skin around the eyes.
3. ageLOC LumiSpa Accent works with specially formulated ageLOC LumiSpa IdealEyes® to speed up and boost benefits to the delicate skin around your eyes, compared to using IdealEyes alone.

What is the rationale behind the name?

“Lumi” stands for luminous and emphasizes the initial benefits

that using the device and products will have on your skin. “Spa” reinforces the fact that the device, treatment heads, and treatment products will deliver spa-quality results every time you use them. “Accent” refers to accentuating your eyes to help them look fresher and brighter.

How does ageLOC LumiSpa Accent work with ageLOC LumiSpa IdealEyes?

ageLOC LumiSpa IdealEyes was specifically developed to use with LumiSpa Accent. They're designed to work together with the oscillating motion of the silicone tip to gently stimulate skin, boosting and increasing the speed of benefits.

Are there any patents on ageLOC LumiSpa Accent?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and ageLOC LumiSpa Accent.

Do I need to use ageLOC LumiSpa IdealEyes with ageLOC LumiSpa, or can I use other eye creams with it?

ageLOC LumiSpa IdealEyes is specially formulated to use with ageLOC LumiSpa. ageLOC LumiSpa IdealEyes enables precise cushioning and interaction with the skin, and helps provide better and faster results compared to using products alone. IdealEyes is the only formulation tested and approved for use with ageLOC LumiSpa Accent to deliver the full system benefits and promote healthy, brighter skin around the eyes.

Do I need to wash my makeup off before using LumiSpa Accent?

Yes. Make sure to use ageLOC LumiSpa Accent after cleansing and toning your face.

Can I use ageLOC LumiSpa Accent on my body?

ageLOC LumiSpa Accent was designed to target the delicate skin around your eyes. We do not recommend you use this attachment anywhere else on your body.

How often should I replace my ageLOC LumiSpa Accent silicone tip?

We recommend you replace your ageLOC LumiSpa Accent tip every three months to ensure you have the best experience.

WHAT'S IN IT

AGELOC LUMISPA IDEALEYES®

Water (Aqua), Glycerin, Dicaprylyl Carbonate, Tribehenin PEG-20 Esters, Dimethicone, Squalane, Cetearyl Alcohol,

ageLOC® LumiSpa® ACCENT

C12-15 Alkyl Benzoate, Cyclopentasiloxane, Octyldodecanol, Butylene Glycol, Tetrahexyldecyl Ascorbate, Caprylic/Capric Triglyceride, Albizia Julibrissin Bark Extract, Caesalpinia Spinosa Fruit Extract, Cyclohexasiloxane, Bisabolol, Geranylgeranylisopropanol, Caffeine, Pentylene Glycol, Tribehenin, Darutoside, Palmitoyl Hexapeptide-12, Palmitoyl Tripeptide-1, Adenosine, Palmitoyl Tetrapeptide-7, Hydrogenated Phosphatidylcholine, PEG-8, Narcissus Tazetta Bulb Extract, Ribes Nigrum (Black Currant) Seed Oil, Carbomer, Cardiospermum Halicacabum Flower/Leaf/Vine Extract, Ceramide NG, Rosmarinus Officinalis (Rosemary) Leaf Extract, Kappaphycus Alvarezii Extract, Helianthus Annuus (Sunflower) Seed Oil, Schizandra Chinensis Fruit Extract, Hydrogenated Lecithin, Sodium Methyl Stearoyl Taurate, Chlorella Vulgaris Extract, Hydroxyphenyl Propamidobenzoic Acid, PEG-10 Phytosterol, Helianthus Annuus (Sunflower) Seed Oil Unsaponifiables, Glyceryl Acrylate/Acrylic Acid Copolymer, Xanthan Gum, Polysorbate 20, Chlorphenesin, Hydroxyacetophenone, Sodium Hydroxide.