

# NU SKIN® FACIAL SPA AND NU SKIN® CONDUCTIVE GEL



## MEET NU SKIN FACIAL SPA AND NU SKIN CONDUCTIVE GEL

Featuring two interchangeable conductors, a large LED display, easy-to-use design, and ergonomic feel, Nu Skin® Facial Spa uses microcurrent technology to stimulate and tone your skin. Indulge in a rewarding 10-minute regimen three times each week for visible results that will make your friends wonder how you find the time to go the spa. **Note:** Please remember to read all directions before using Facial Spa.

### WHY YOU'LL LOVE IT

- Stimulates skin for improved appearance.
- Provides facial toning—improving facial contour and making skin look more voluminous.
- Features microcurrent technology.
- Includes two interchangeable conductors that are specially designed to follow the contours of the face while keeping the Conductive Gel in constant contact with your skin.
  - Large Conductor—designed for use during full-face or forehead regimens.
  - Small Conductor—helps deliver current to specific areas during targeted area regimens.
- Designed with a modern, sleek user interface.

### WHAT IT WORKS WITH

Nu Skin Conductive Gel is specifically formulated for use with Nu Skin Facial Spa to ensure adequate conductivity throughout your spa session and help Nu Skin Facial Spa glide across skin.

**Ingredients:** Water (Aqua), Glycerin, Pentylene Glycol, Carbomer, Sodium Hydroxide, Chlorphenesin.

## LEARN MORE ABOUT IT

### Is Nu Skin Facial Spa a medical device?

Nu Skin has received FDA clearance to market Nu Skin Facial Spa as a class II medical device for over-the-counter cosmetic use. Facial Spa with its accompanying Conductive Gel is a facial toning device designed to stimulate and tone the skin for an improved appearance.

### What are the session times?

Nu Skin Facial Spa features a 10-minute, full-face, basic regimen, and three five-minute area regimens. The basic regimen has a five-minute Spa Session A and a five-minute Spa Session B. All of the area regimens include a two-minute Spa Session A and a three-minute Spa Session B.

### How often should I use Nu Skin Facial Spa?

You can conduct a Nu Skin Facial Spa regimen up to three times each week. Each regimen consists of two sessions—Session A and Session B. Facial Spa should be only used as directed with Conductive Gel.

### Can Facial Spa be used on the neck?

Nu Skin Facial Spa is registered for facial use only and should not be used on other parts of the body, including the neck.

### Why does Facial Spa beep?

When you begin a session, the device display will be lit and the device will beep every 10 seconds to indicate that the unit is on and working. The audible beeps are the **ONLY** indicator that the device is working. Once the spa session time reaches 30 seconds remaining, the device will beep every five seconds until the spa session is complete, at which point the “spa session complete” alert will sound.

# NU SKIN® FACIAL SPA AND NU SKIN® CONDUCTIVE GEL

## How long will the batteries last?

In-house testing showed that batteries lasted several months with normal usage.

## WARNINGS

To reduce the risk to your health: use product only for its intended use, as described in the Facial Spa User Manual.

## CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND ADVERSE REACTIONS

⊗ Please read the following advisories before using your Facial Spa

### CONTRAINDICATIONS

Do not use the Facial Spa if you:

- Have a metal implant or electronic implanted device.
- Have suspected or diagnosed epilepsy or suffer from seizures.
- Have metal braces on your teeth.
- Are allergic to gels or metal or are sensitive to contact with chrome.
- Are pregnant or attempting to become pregnant.
- Have open sores or wounds.
- Are prone to rosacea.
- Have problematic or sensitive skin.

The Facial Spa device is not intended for use by persons (including children) with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge.

The Facial Spa and Nu Skin Conductive Gel are contraindicated for use in the following areas:

- Eyes, eyelids, and surrounding tissue (orbicularis oculi muscle);
- Lips (labium superius oris and labium inferius oris);
- Throat and chest; and
- Any areas that cannot sense touch.

### WARNINGS

- ⊗ If you are in the care of a physician, consult with your physician before using the device.
- ⊗ Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- ⊗ Do not apply stimulation across your chest because the introduction of electrical current into the chest

may cause rhythm disturbances to your heart, which could be lethal.

- ⊗ Do not apply stimulation over painful areas. If you have painful areas, you should consult with your physician before using the device.
- ⊗ Do not apply stimulation over open wounds or rashes; swollen, red, infected, or inflamed areas; or over skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- ⊗ Do not apply stimulation over, or in proximity to, cancerous lesions or diseased tissue.
- ⊗ Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms) which may not operate properly when an electrical stimulation device is in use.
- ⊗ Do not apply stimulation while sleeping.
- ⊗ Do not change, release, or attach conductors when the device is turned ON. You risk exposing yourself to the exposed electrical lead, which may result in electrical shock or burnt skin.
- ⊗ Do not conduct more than three Facial Spa regimens per week.
- ⊗ Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- ⊗ Adult use only—do not use the device on children; it has not been evaluated for pediatric use.
- ⊗ Consult with your physician before using the device. The device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- ⊗ Apply stimulation only to normal, clean, healthy skin.
- ⊗ Do not operate while bathing or showering. To reduce the risk of shock, burn, fire, or injury, do not immerse your device in water. Do not place or store the device where it may fall or be pulled into a bathtub, shower, sink, or toilet.
- ⊗ Do not apply stimulation on skin without Nu Skin Conductive Gel. Using the device without Conductive Gel may burn your skin. Only use with Nu Skin Conductive Gel.

### PRECAUTIONS

- ⊗ Use the Facial Spa device only as directed.
- ⊗ Always inspect your Facial Spa and conductors before use. Ensure there are no rough or defective surfaces, sharp edges, blistering of the chrome, and/or abnormalities on any surface. If you notice any abnormalities, do not use the device and contact Nu Skin Support Services.

# NU SKIN® FACIAL SPA AND NU SKIN® CONDUCTIVE GEL

- ⊘ The long-term effects of electrical stimulation are unknown.
- ⊘ Since the effects of stimulation on the brain are unknown, stimulation should not be applied across your head, and electrodes should not be placed on opposite sides of your head.
- ⊘ The safety of electrical stimulation during pregnancy has not been established.
- ⊘ You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (Nu Skin Conductive Gel).
- ⊘ If you have suspected or been diagnosed with heart disease, you should follow precautions recommended by your physician.
- ⊘ Use caution if you have a tendency to bleed internally, such as following an injury or fracture.
- ⊘ Consult with your physician prior to using the device after a recent surgical procedure because stimulation may disrupt the healing process.
- ⊘ Use caution if stimulation is applied over the menstruating uterus.
- ⊘ Do not use over areas of skin that lack normal sensation, as it could result in burns.
- ⊘ Keep the device out of the reach of children. Children should be supervised to ensure that they do not play with the device.
- ⊘ Use the device only with the conductors and accessories recommended by Nu Skin.
- ⊘ Always turn device off (display will appear unlit) before removing/replacing conductor and do not use without either Large Conductor or Small Conductor locked into place; otherwise there is risk of shock.
- ⊘ Never operate the device if it has been damaged. Using a damaged device may put you at risk of shock, burn, fire, or injury; or may permanently damage the device.
- ⊘ Do not apply excessive pressure during use. Only apply enough gentle pressure to maintain constant contact between the conductor and your skin.
- ⊘ Users may experience a mild oral sensation during use. Discontinue use if you experience any discomfort.
- ⊘ Discontinue use and call Nu Skin Support Services if you experience stinging, redness, or irritation.
- ⊘ Do not use for longer periods of time or more frequently than recommended.

## ADVERSE REACTIONS

- You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin.
- You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face. Discontinue use if this occurs.

Discontinue use and report any adverse reactions to Nu Skin Support Services. If you have any questions related to the safe use of this device, please call Nu Skin Support Services.

Nu Skin Support Services  
1-800-487-1000